



Power Range Guidelines

Calgary HD, Calgary World Wide & Digital Superior
SV & Occupational

	All Progressive designs	Standard	Wrap Compensated SV	Occupational
Materials	all Corridors	Aspheric SV	4.00-8.00 Base Only	Office Lens
1,5	+9.00 to -8.50	+9.00 to -12.00	+7.00 to -6.00	+5.00 to -6.00
Trivex	+8.00 to -8.50	+8.00 to -10.00	+7.00 to -6.00	+5.00 to -6.00
Mid-Index	+9.00 to -9.00	+9.00 to -12.00	+7.00 to -6.00	+5.00 to -6.00
Poly	+9.00 to -8.50	+9.00 to -14.00	+7.00 to -6.00	+7.00 to -9.00
1,6	+9.00 to -10.00	+9.00 to -14.00	+7.00 to -6.00	+8.00 to -6.00
1,67	+12.00 to -11.00	+12.00 to -16.00	+7.00 to -6.00	+8.00 to -6.00
1,74	+12.00 to -12.00	+12.00 to -20.00	+7.00 to -6.00	+8.00 to -6.00

Combined Powers

These ranges represent combined power limits - see examples below:

Positive Cylinders in the Distance region = Sphere + Cylinder
 Positive Cylinders in the Near region = Sphere + Cylinder + Add
 Negative Cylinders in the Distance region = Sphere
 Negative Cylinders in Near region = Sphere + Add

Sphere	Cylinder	Add	Combined Power	
			Distance	Near
-2.00	-1.00	2.0	-2.00	Zero
+2.00	+1.00	1.0	+3.00	+4.00
-2.00	+1.00	3.0	-1.00	+2.00
+3.00	-2.00	4.0	+3.00	+7.00

Adds Adds are currently available from +1.00 to +4.50 *These ranges are general guidelines.

Cylinders Cylinders -.12 to -5.00 all products *These ranges are general guidelines.

Materials Ranges will be reduced by all Polarized, Transitions and Photochromics lenses.

Calgary HD, Calgary World Wide and Digital Superior designs have no range limits but the general limit is around -14D although it depends on each specific device.

Sphere	Cylinder	Add	Total	
-11.00	+3.00	0.00	-11.00	OK
-11.00	-3.00	0.00	-14.00	At limit
-11.00	-3.00	2.00	-12.00	OK

Base Curves & Blank Selection	On positive (+) lenses, the bioconvexity limit can be upto 10 Base Curve. But the Sphere and the Add must total less than the Base Curve. E.G. Sphere +4 with +3 Add < Base curve
--	--

